

20210508

Notes on multipoint Z series.

To set up multiple positions (XY) with different Z positions, use the Include Z checkbox in the ND Acquisition window. Use PFS is optional.

(For objects close to the coverslip, such as cultured cells, where asymmetric Z is used to go in the up direction only, PFS may be used instead of Z.)

Z Device in XY tab always set for Ti2 Drive. This sets the Z position absolutely using the motor inside the microscope stand.

In the Z window, you must use either the symmetric or asymmetric mode defined by range. The Relative button **should be grayed out** and **DO NOT click the Home button**.

The reference Z position is defined in the XY tab.

Regardless whether you use the Ti2 Drive, NIDAQ Piezo Z, or Triggered NIDAQ Piezo Z, these operate relative to the absolute Z position set in the XY tab.


When you are setting the Z position, Piezo Z in lower right should always stay at 110 um.

Devices > Move Piezo Z to Home Position

will set to 110 um and

Devices > Keeps Z position and centers Piezo Z

should lock it at 110 um except during Piezo Z specific addressing.



Always click Split Multipoints at bottom of the Advanced menu in XY tab so that each position will be saved separately. At the end of the experiment, click No to the question about splitting multipoints as they are already saved separately on disk. Doing this 1.) allows you to access in ImageJ the data during image collection, 2.) makes smaller files, and 3.) saves time at the end of the experiment because you do not need to split the multipoints.

ND Acquisition x

Experiment: ND Acquisition

M:

Z:

Save to File

Path:

Filename:

Custom Metadata

Order of Experiment:

Time XY Z Large Image λ

Points Move Stage to Selected Point

Point Name	X [mm]	Y [mm]	Z [μm]	
<input checked="" type="checkbox"/> #1	-6.426	0.694	3102.060	
<input checked="" type="checkbox"/> #2	-> -6.711	0.736	3107.460	<- Offset All X,Y,Z
<input checked="" type="checkbox"/> #3	-5.996	0.625	3108.060	
<input type="checkbox"/>				

Include Z Relative XY

Close Active Shutter during Stage Movement Use PFS

Z Device:

Redefine Reference Z after Auto Focus

Leave PFS offset ON between points

Autofocus: Use Focus Surface

Execute Command before Capture

Execute Command after Capture

Split Multipoints

ND Acquisition x

Experiment: ND Acquisition

M:

Z:

Save to File

Path: E:\MCammer E

Filename: 202100408_test003.nd2

Custom Metadata

Order of Experiment


Time XY Z Large Image λ

Top

Relative Home

Bottom

symmetric



Step: 6 μm 0.4 μm 7 Steps Range: 36.00 μm

Range: 36 μm <-18.00, +18.00> Relative Positions: Top: +18.00 μm

Z Device: Triggered NIDAQ Piezo Z Piezo Bottom: -18.00 μm

Close Active Shutter during Z Movement Direction: Bottom to Top Top to Bottom

Execute Command before Capture

Execute Command after Capture

ND Acquisition x

Experiment: ND Acquisition

M:

Z:

Save to File


Path: E:\MCammer E

Filename: 202100408_test003.nd2

Custom Metadata

Order of Experiment

Time XY Z Large Image λ



Step: 6 μm \leftarrow 0.4 μm 7 Steps Range: 36.00 μm

Below: -0.00 μm Above: +36.00 μm Relative Positions:

Z Device: Triggered NIDAQ Piezo Z Top: +36.00 μm

Bottom: -0.00 μm

Close Active Shutter during Z Movement Direction: Bottom to Top Top to Bottom

Execute Command before Capture

Execute Command after Capture